BLUE LION
LIFESTYLE

BLUE LION
LIFESTYLE

KOSEY

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Congratulations!

You must be proud that you have invested in yourself and your physical health. When you are working hard every day, you also need to take good care of your body. In this user guide you can read how you can safely give your back a nice and loose massage by your new personal back masseur, so that you can get after it again!

About your new personal back masseur: Kosey

Your new back massager, called Kosey, solves a problem that many people face these days. More and more people are experiencing problems with their backs and sitting a lot does not help. Until now. It is Kosey's mission to help you improve your physical health. Put him in your favorite seat and experience a deep massage from the four big massage heads that loosen up your back muscles in three different zones:

1. Full back
2. Upper back
3. Lower back

Kosey massages you in the zone you have selected. Because many people mainly suffer from pain in their lower back and around their shoulderblades, he massages extra deeply in these areas. Do you have a specific spot where you want to be massaged? If you press the 'spot' button on the remote control, he will focus his massage on a specific place and you can manually guide him up or down with the arrows. You can determine the massage intensity with the remote control and by using your massage cushion with or without a protective cover. You can also put an extra towel over it or put on a thick sweater if you want to make it even less intense. The seat has a vibration option with three different intensities so you can enjoy those 'good vibes'.

Use your back massage cushion on the couch, in the chair or in the car. Loosen your muscles before a workout, promote recovery after an injury, stimulate your circulation while working in the office, or enjoy a relaxing massage after a long day at work. However you choose to use it, we hope you're as excited about it as we are!

Before you start...

Read the safety instructions on the next page. We understand that you are excited and want to use your massage cushion right away, so we have kept it as short as possible to save your time and also save paper.

Warning

Careless use of the massage cushion can potentially cause personal injury and damage to the device or the environment. Read the safety instructions carefully to continue to enjoy your back massage cushion in the long term.

Safety instructions

1. This device is not suitable for animals or use while driving (only use in a stationary car)
2. Only use your your back massager after consulting a doctor or specialist if:
 - you are pregnant;
 - you have a pacemaker;
 - you have artificial joints or implants;
 - you experience intense pain during the massage or have physical health conditions such as: heat sensitivity, blood circulation problems, varicose veins, bruises, muscle or tissue inflammation, injuries and open wounds.
3. Unplug the massager after every use.
4. Grab and carry your back massager by itself: do not lift or pull on the cable.
5. Avoid your back massager from getting wet and only handle it with dry hands.
6. Only use the original adapters and power cable supplied by the manufacturer.
7. Keep your back massager and its power supplies away from heat sources, highly inflammable substances, vapour, liquids or sharp objects.
8. Before use, check if the cables are free from knots, twists or damage.
9. Provide adult supervision when a child, elderly or disabled person uses the back massager.
10. Clean with a soft absorbent cloth without detergents (see p.8).
11. For your own safety, parts should only be replaced by professional / authorized service personnel.
12. Send the device to a professional service center for repair and maintenance.
13. Lay the cable carefully so that no one can trip over it.
14. Your back massager is intended for indoor use.
15. Your back massager is intended for private use.
16. Keep all parts and packaging materials out of the reach of children.

Important to know:

The warranty is void if these safety instructions are not observed or if the back massager is used for purposes other than those described above. The supplier cannot be held liable for damage in these cases.

What's in the box?

1x Blue Lion Kosey-1 Shiatsu back massager

1x Remote control

1x 12 V Home adapter

1x 12 V Car adapter

1x User guide

Technical details

Brand	Blue Lion
Model	Kosey-1
Voltage Supply	100-240 V - 50/60 Hz
Power Supply Unit	12 V 1,5 A
Power Consumption	18 W
Timer	15 minutes
Dimensions	40.9 x 18.1 inch
Weight	8.6 lbs
Power cable length	110.2 inch
Car adapter cable length	70.9 inch
Test reports and certifications	CE, RoHS, LVD, EMC, ERP

RoHS

Did you know...

The massage cushion is named after a mix of the word 'cozy' and an African name 'Kosey', which means: Lion


8 Steps to your perfect massage experience

- 1 Place your Blue Lion back massager on your favorite chair or sofa.
- 2 Use the elastic straps on the back of your massager to attach it to the chair for extra support.
- 3 Connect the adapter cable to the cable attached to the cushion as shown in the picture below:

- 4 Plug the home adapter into an electric outlet.
- 5 Turn on your back massager with the remote control and choose your favorite zone and mode (see controller instructions on next page).
- 6 Choose your preferred massage intensity.
 - Too intense? Use the detachable protection cover for a gentler massage.
 - Not intense enough? Remove the detachable protection cover for a more intense massage.

Is the massage at the lowest intensity with the protection cover still too intense? You can use an extra towel or blanket or put on a nice thick hoodie to make it more soothing.

- 7 If necessary, use a seat cushion to adjust the height of your back.
- 8 Put the remote control in the pouch on the right side of the seat and enjoy your personalized massage!

Controller instructions

Massage zones

Where do you want to be massaged?

- Full back
- Upper back
- Lower back

Select your favorite massage zone.

Vibration


Activate vibration in the seat:

- L= low
- M= Medium
- H= High

Intensity

Select your favorite massage-intensity:

- L= Low
- H= High

Spot massage

Press 'Spot' to keep Kosey massage that specific spot. Use the arrows to manually guide him up or down.

Heat

Activate the heating in the massage heads. (perfect for after a hard workout).

Demo

Quick demonstration of all different functions. The flashing lights indicate the active functions. The back massager automatically turns off after completion.

Power

Switch your back-massager on or off.

Kosey will automatically turn off after 15 minutes.

Do you want Kosey to continue massaging? Press the power button again.

Cleaning

Take good care of Kosey and he will take good care of you. Careless maintenance can shorten its life and void the warranty. This is how to keep Kosey clean:

1. Unplug the massage cushion.
2. Use a soft absorbent (slightly damp) cloth.
3. Immediately dry your back-massager cushion afterwards.
4. Only use your massager again until he is completely dry.

✓ Service & Warranty

We believe in the quality of our Kosey-1 back-massager with great confidence, so we offer you a two year warranty against defects due to material or manufacturing defects. Please keep your proof of purchase and this manual safe in case you need it.

Note: The warranty only applies if the user has carefully followed the instructions in this manual. The warranty does not apply to damage caused by improper/careless use, or wear and tear. Read the extended warranty conditions on our website for more information.

Contact

Do you have any questions, comments or just want to share something with us? You can call or mail us and send us a WhatsApp message too. Whatever it is, we are here for you.

help@bluelionlifestyle.nl

085-0604658

www.bluelionlifestyle.nl

Help nature grow

We hope you enjoy Kosey's massages for many years to come. Has your back-massager reached the end of its life after a long enjoyable time?

Don't just take it to the dump, but hand it in at a recycle collection point and be a hero! By recycling the cardboard, plastic and equipment, together we can contribute to the growth of nature and a cleaner living environment.

Share your photos, video's and stories about your back-massager with us!

Tag us @bluelionlifestyle and we'll share your post with our followers so you can help and inspire each other.

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© 2026 Blue Lion Lifestyle